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UCL

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BD Phoenix™ ID/AST Manual Panel Inoculation

McFarland 0.5 ~0.6

- 1** Open panel package.
INOCULATE WITHIN
2 HOURS.

- 4** After panel inoculation, place closures
securely on the panel to seal.

Prepare Purity Plate.
Using a sterile loop,
recover a small drop from
the inoculum fluid either
before or after inoculating
the panel. Inoculate an agar
plate (any appropriate medium)
for purity check. Incubate plates
for 18 - 24 hours or 18 - 48 hours
for yeast, at 35°C under appropriate
conditions.

- 2** Prepare Phoenix™ ID Broth.
Inoculate panel with prepared
ID Broth within (≤) 60 MINUTES.

PREPARE ID BROTH.

- Add organism to ID Broth.
- Vortex.
- Make .50 - .60 McFarland (Standard inoculum)*.
- OR
- Make .20 - .30 McFarland (Low inoculum)**.
- OR
- Make 2.00 - 2.40 McFarland (Yeast inoculum).

- 5** Load panel(s) into BD Phoenix instrument
within (≤) 30 MINUTES after inoculation.

- 3** Prepare Phoenix™ AST or AST-S Broth.
(ID/AST, AST Panels only)
Inoculate panel with prepared broth within (≤) 30 MINUTES.

PREPARE AST or AST-S BROTH***

- Add one free falling drop of AST or AST-S Indicator**** into AST or AST-S Broth. Tube can be used for up to 2 hours if exposed to light or up to 8 hours in the dark. Mix by inverting.
- Add 25 µL prepared ID Broth (Standard inoculum).
- OR
- Add 50 µL prepared ID Broth (Low inoculum).

BD Diagnostics
7 Loveton Circle
Sparks, MD 21152-0999
800.638.8663
www.bd.com/us

BD, BD logo and Phoenix are trademarks of
Becton, Dickinson and Company. ©2012 BD
2-222341 January 2012
SHP-228

BD Phoenix™ ID/AST

* Not applicable for yeast ID panels.
** Not applicable for strep or yeast ID panels.
*** For gram-negative or gram-positive panels
use Phoenix AST Indicator and AST Broth.
For strep panels use Phoenix AST-S Indicator and AST-S Broth.
**** AST or AST-S indicator is good for 14 days after opening (refrigerated).